
UofA Laboratory Documentation

Release 0.4

Sam Jennings

Feb 18, 2019

CONTENTS:

| | | |
|----------|--------------------------------|-----------|
| 1 | Package Laboratory | 1 |
| 1.1 | Module Config | 1 |
| 1.2 | Module Setup | 1 |
| 2 | Package Utils | 5 |
| 2.1 | Module Calibrate | 5 |
| 2.2 | Module Data | 5 |
| 2.3 | Module Loggers | 7 |
| 2.4 | Module Notifications | 7 |
| 2.5 | Module Plotting | 7 |
| 3 | Package Drivers | 9 |
| 3.1 | Module DAQ | 9 |
| 3.2 | Module Furnace | 10 |
| 3.3 | Module LCR | 12 |
| 3.4 | Module MFC | 12 |
| 3.5 | Module Motor | 15 |
| 3.6 | Module Other | 16 |
| 4 | Indices and tables | 17 |
| | Python Module Index | 19 |
| | Index | 21 |

PACKAGE LABORATORY

1.1 Module Config

This is a configuration file for setting up the laboratory. It contains settings to set the name of the the experiment, sample dimensions, instrument addresses and physical constants. Experiment name and sample dimensions should be modified with each new experiment. Everything else can remain as is unless the physical setup of the lab has changed.

1.2 Module Setup

```
class laboratory.setup.Setup (filename=None, debug=False)
```

Bases: object

Sets up the laboratory

| Attributes | Description |
|------------|--|
| data | houses lab data during measurements or after parseing |
| logger | creates a logger for error reporting |
| plot | contains different plotting tools for data visualisation |
| command | contains the drivers for controlling instrumentation |

| Methods | Description |
|------------------|--|
| device_status | checks the status of all connected devices |
| get_gas | retrieves and saves data from a single mass flow controller |
| get_impedance | retrieves and saves complex impedance data from the LCR meter |
| get_temp | retrieves and saves temperature data from the furnace |
| get_thermo | retrieves and saves thermopower data from the DAQ |
| load_data | will parse a datafile and store in 'data' structure for post- processing and visualisation |
| load_frequencies | loads a set of frequencies into Data object |
| load_instruments | connects to all available instruments |
| run | begins a new set of laboratory measurements |

Example

```
>>> import Laboratory
>>> lab = Laboratory.Setup()
>>> lab.run('some_controlfile')
```

append_data (filename)

debug

delayed_start

Starts the experiment at a given time the next day. Can be set to any 24 hour time in string format.

Example

```
>>> import Laboratory
>>> lab = Laboratory.Setup()
>>> lab.delayed_start = '0900' #start at 9am the following day
>>> lab.run('some_controlfile')
```

device_status()

Checks the status of all devices. If desired, this function can send an email when something has become disconnected

Returns True if all devices are connected and False if any are disconnected

Return type Boolean

get_gas(gas_type)

Gets data from the mass flow controller specified by gas_type and saves to Data structure and file

Parameters **gas_type** (*str*) – type of gas to use when calculating ratio (either 'h2' or 'co')

Returns [mass_flow, pressure, temperature, volumetric_flow, setpoint]

Return type list

get_impedance()

Sets up the lcr meter and retrieves complex impedance data at all frequencies specified by Data.freq. Data is saved in Data.imp.z and Data.imp.theta as a list of length Data.freq. Values are also saved to the data file.

get_temp(target)

Retrieves the indicated temperature of the furnace and saves to Data structure and file

Note: this is the temperature indicated by the furnace, not the temperature of the sample

Parameters **target** (*float*) – target temperature of current step

get_thermopower()

Retrieves thermopower data from the DAQ and saves to Data structure and file

Returns [thermistor, te1, te2, voltage]

Return type list

load_data(filename)

loads a previous data file for processing and analysis

Parameters **filename** (*str*) – path to data file

load_frequencies(min=20, max=2000000, n=50, log=True, filename=None)

Loads an np.array of frequency values specified by either min, max and n or a file containing a list of frequencies specified by filename.

Parameters

- **filename** (*str*) – name of file containing frequencies
- **n** (*int*) – number of desired frequencies

- **min** (*int, float*) – minimum frequency (Hz) - may not be below default value of 20 Hz
- **max** (*int, float*) – maximum frequency (Hz) - may not exceed default value of 2×10^6 Hz
- **log** (*boolean*) – specifies whether array is created in linear or log space. default to logspace

Example

```
>>> lab = Laboratory.Setup()
>>> lab.load_frequencies(min=1000, max=10000, n=10)
>>> print(lab.data.freq)
[1000 2000 3000 4000 5000 6000 7000 8000 9000 10000]
>>> lab.load_frequencies(min=1000, max=10000, n=10, log=True)
>>> print(lab.data.freq)
[1000 1291.55 1668.1 2154.43 2782.56 3593.81 4641.59 5994.84 7742.64 10000]
```

load_instruments()

Loads all the laboratory instruments. Called automatically when calling Setup() without a filename specified.

Returns lcr, daq, mfc, furnace, motor

Return type instrument objects

preflight_checklist (*controlfile*)

Conducts necessary checks before running an experiment abs

Parameters **controlfile** (*string*) – name of control file for the experiment

reconnect()

Attempts to reconnect to any instruments that have been disconnected

restart_from_backup()

TODO - reload an aborted experiment and pick up where it left off

run (*controlfile=False*)

starts a new set of measurements. requires a control file that contains specific instruction for the instruments to follow. see the tutorial section for help setting up a control file.

Parameters **controlfile** (*str*) – path to control file

save_data (*val_type, vals, gastype=None*)

Takes input values and saves to both the current Data object and an external file

Parameters

- **val_type** (*str*) – type of measurement being saved
- **vals** (*list*) – the required values suitable to that specified by val_type
- **gastype** (*str*) – [optional] required when saving gas data

set_fugacity (*buffer, offset, gas_type*)

Sets the correct gas ratio for the given buffer. Percentage offset from a given buffer can be specified by 'offset'. Type of gas to be used for calculations is specified by gas_type.

Parameters

- **buffer** (*str*) – buffer type (see table for input options)
- **offset** (*float, int*) – percentage offset from specified buffer

- **gas_type** – gas type to use for calculating ratio - can be either 'h2' or 'co'

shut_down()

Returns the furnace to a safe temperature and closes ports to both the DAQ and LCR. (TODO need to close ports to motor and furnace)

PACKAGE UTILS

2.1 Module Calibrate

`laboratory.utils.calibrate.find_center(self)`

TODO - Attempts to place the sample at the center of the heat source such that $te1 = te2$. untested.

`laboratory.utils.calibrate.furnace_profile()`

Records the temperature of both electrodes ($te1$ and $te2$) as the sample is moved from one end of the stage to the other. Used to find the center of the stage or the $xpos$ of a desired temperature gradient when taking thermopower measurements.

2.2 Module Data

class `laboratory.utils.data.Data` (*freq=None, filename=None*)

Bases: `object`

Storage for all data collected during experiments. Data file are loaded into this object for processing and plotting

| Attributes | Description |
|-----------------------|---|
| <code>freq</code> | array of frequencies for use by the LCR meter |
| <code>filename</code> | name of the file being used |
| <code>time</code> | times for each measurement |
| <code>thermo</code> | stores thermopower data |
| <code>gas</code> | stores gas and fugacity data |
| <code>temp</code> | stores temperature data |
| <code>imp</code> | stores impedance data |
| <code>xpos</code> | stores stage x position at each measurement |

Example

```
>>> lab = Laboratory.Setup('somefile.dat')
>>> print(lab.data.temp.indicated)
[100, 105, 110, 115, 120]
```

filename

freq

class `laboratory.utils.data.Gas`

Bases: `object`

Stores the separate gas data under one roof

| Attributes | Description |
|------------|----------------------------------|
| h2 | hydrogen flow rate |
| co2 | carbon dioxide flow rate |
| co_a | carbon monoxide corase flow rate |
| co_b | carbon monoxide corase flow rate |

class `laboratory.utils.data.Impedance`

Bases: `object`

Stores complex impedance data

| Attributes | Description |
|------------|-------------|
| Z | impedance |
| theta | phase angle |

class `laboratory.utils.data.MFC_data`

Bases: `object`

Stores gas data for an individual mass flow controller

class `laboratory.utils.data.Temp`

Bases: `object`

Stores furnace temperature data

| Attributes | Description |
|------------|-------------------------------------|
| target | target temperature of current cycle |
| indicated | temperature indicated by furnace |

class `laboratory.utils.data.Thermo`

Bases: `object`

Stores thermopower data

| Attributes | Description |
|------------|--|
| tref | temperature of the internal thermistor |
| te1 | temperature of electrode 1 |
| te2 | temperature of electrode 2 |
| volt | voltage across the sample |

`laboratory.utils.data.append_data(filename, data)`

`laboratory.utils.data.load_data(filename)`

`laboratory.utils.data.save_obj(obj, filename)`

Saves an object instance as a .pkl file for later retrieval. Can be loaded again using `:meth:'Utils.load_obj'`

Parameters

- **obj** (*class*) – the object instance to be saved
- **filename** (*str*) – name of file

2.3 Module Loggers

`laboratory.utils.loggers.data()`

Sets up the data file in much the same way as the log file. Data cannot be output to the console. Data file can be found in /datafiles/

`laboratory.utils.loggers.lab(name)`

Sets up logging messages for the laboratory. Sends to both a file and the console by default. Levels for both the file and console can be set to anything defined by the python logging package (DEBUG, INFO, WARNING, ERROR, CRITICAL). Specified log level AND GREATER will be included. Logfiles can be found in /logfiles/

2.4 Module Notifications

`class laboratory.utils.notifications.Messages`

Bases: object

`delayed_start = "It just ticked over to {} so I'm going to set up the instruments and"`

`device_error = "I've got some bad news! The {} is no longer sending or receiving messa"`

`step_complete = "Just letting you know that step {} is now complete! I'm going to set t"`

`laboratory.utils.notifications.send_email(toaddr, message, cc=False, logfile=False, datafile=False)`

Sends an email to the specified email address. logfile or datafile can be attached if desired. used mainly for email updates on progress during long measurement cycles. mailer is geophysicslabnotifications@gmail.com.

Parameters

- **toaddr** (*str*) – full email address of intended recipient
- **message** (*str*) – message to include in email
- **cc** (*str, list*) – email can be carbon copied to additional addresses in cc
- **logfile** (*boolean*) – whether to attach the current logfile
- **datafile** (*boolean*) – whether to attach the current datafile

2.5 Module Plotting

Contains plotting tools for use with laboratory data.

| Class Objects | Description |
|---------------|--|
| LabPlots | houses an assortment of plotting tools |

| Methods | Description |
|------------------------------------|---|
| <code>impedance_fit</code> | estimates the diameter of the impedance arc based on |
| <code>dt_to_hours</code> | converts recorded datetimes to time elapsed |
| <code>leastsq_circle</code> | fits a least squares circle to impedance data |
| <code>index_temp</code> | index the nearest temperature value |
| <code>get_Re_Im</code> | returns the real and imaginary components of complex impedance |
| <code>calculate_resistivity</code> | calculates resistivity from impedance spectra and sample dimensions |

```
class laboratory.utils.plotting.LabPlots(data)
```

Bases: object

Contains an assortment of useful plots for visualising with the laboratory data object

| Attributes | Description |
|---------------------------|---|
| <code>data</code> | the data object from the laboratory |
| <code>time_elapsed</code> | time elapsed since the start of the experiment [in hours] |

| Methods | Description |
|----------------------------|--|
| <code>arrhenius</code> | plots log conductivity vs reciprocal temperature [K] |
| <code>cond_time</code> | plots conductivity vs elapsed time |
| <code>cole</code> | creates a cole-cole plot at a given temperature/s |
| <code>gas</code> | plots massflow vs elapsed time |
| <code>temperature</code> | plots tref, te1, te2, and target vs elapsed time |
| <code>imp_diameter</code> | plots imp_diameter vs elapsed time |
| <code>cond_fugacity</code> | plots log conductivity vs fugacity |
| <code>voltage</code> | plots voltage vs elapsed time |

```
arrhenius ()
```

Plots inverse temperature versus conductivity

```
cole (temp_list, start=0, end=None, fit=False)
```

Creates a Cole-Cole plot (imaginary versus real impedance) at a given temperature. Finds the available self.data to the temperature specified by 'temp'. A linear least squares circle fit can be added by setting fit=True.

Parameters `temp` (*float/int*) – temperature in degrees C

```
cond_fugacity ()
```

Plots inverse temperature versus conductivity

```
cond_time ()
```

Plots conductivity versus time

```
gas ()
```

Plots mass_flow self.data for all gases versus time elapsed

```
imp_diameter ()
```

Plots the impedance diameter versus time_elapsed

```
temperature ()
```

Plots furnace indicated and target temperature, thermocouple temperature and thermistor self.data versus time elapsed

```
voltage ()
```

Plots voltage versus time

```
laboratory.utils.plotting.calculate_resistivity (Z, theta)
```

Calculates the resistivity of the sample from the resistance and sample dimensions supplied in config.py

```
laboratory.utils.plotting.dt_to_hours (start, dlist)
```

```
laboratory.utils.plotting.get_Re_Im (Z, theta)
```

```
laboratory.utils.plotting.impedance_fit (Z, theta)
```

```
laboratory.utils.plotting.index_temp (T, Tx)
```

```
laboratory.utils.plotting.leastsq_circle (x, y)
```

PACKAGE DRIVERS

3.1 Module DAQ

class laboratory.drivers.daq.DAQ

Bases: object

Driver for the 34970A Data Acquisition / Data Logger Switch Unit

| Attributes | message |
|------------|--|
| maxtry | max number to attempt command |
| status | whether the instrument is connected |
| therm | specifies type of thermistor |
| tref | '101' - channel for thermistor |
| te1 | '104' - channel for electrode 1 |
| te2 | '105' - channel for electrode 2 |
| volt | '103' - channel for voltage measurements |
| switch | '205', '206' - channels for switch between LCR and temp measurements |
| address | computer port address |

| Methods | message |
|---------------|---|
| connect | attempt to connect to the LCR meter |
| configure | configures device for measurements |
| get_temp | gets temperature from te1,te2 and tref |
| get_voltage | gets voltage measurement |
| read_errors | reads errors stored in the DAQ |
| reset | resets the device |
| shut_down | shuts down the device |
| toggle_switch | switches configuration between temp and voltage |

Note: do not change class attributes unless the physical wiring has been changed within the DAQ

configure()

Configures the DAQ according to the current wiring

connect()

Connects to the DAQ

get_temp()

Scans the thermistor and thermocouples for temperature readings

Returns [tref,te1,te2]

Return type list of floats (degrees Celsius)

get_thermopower ()

Collects both temperature and voltage data and returns a list

get_voltage ()

Gets voltage across the sample from the DAQ

Returns voltage

Return type float

read_errors ()

Reads errors from the DAQ (unsure if working or not)

reset ()

Resets the device

shutdown ()

Shuts down the DAQ

toggle_switch (*command*)

Opens or closes the switch to the lcr. Must be closed for impedance measurements and open for thermopower measurements.

Parameters **command** (*str*) – either ‘thermo’ to make thermopower measurements or ‘impedance’ for impedance measurements

3.2 Module Furnace

class `laboratory.drivers.furnace.Furnace` (*ports=None*)

Bases: `object`

Driver for the Eurotherm 3216 Temperature Controller

Note: units are in °C

| Attributes | message |
|--------------|---|
| maxtry | max number to attempt command |
| default_temp | revert to this temperature when resetting |
| status | whether the instrument is connected |
| address | computer port address |

| Methods | message |
|--------------|-------------------------------------|
| connect | attempt to connect to the LCR meter |
| set_temp | set temperature of furnace |
| get_temp | get temperature from furnace |
| set_heatrate | set heatrate of furnace |
| get_heatrate | get heatrate from furnace |
| set_other | set another parameter on furnace |
| get_other | get another parameter from furnace |
| reset | resets the device |

configure()

flush_input()

flush_output()

heating_rate (*heat_rate=None, address=35*)
 Sets the desired heating rate of furnace. Modbus address - 35

Parameters **heat_rate** (*float, int*) – heating rate in °C/min

Returns True if succesful, False if not

Return type Boolean

indicated (*address=1*)
 Query current temperature of furnace. Modbus address - 1

Returns Temperature in °C if succesful, else False

Return type float/boolean

other (*address, value=None*)
 set value at specified modbus address.

Parameters

- **modbus_address** (*float, int*) – see furnace manual for addresses
- **val** (*float, int*) – value to be sent to the furnace

Returns True if succesful, False if not

Return type Boolean

reset()
 resets the furnace to default temperature

setpoint_1 (*temperature=None, address=24*)
 Sets target temperature of furnace. Modbus address - 24

Parameters **temp** (*float, int*) – temperature in °C

Returns True if succesful, False if not

Return type Boolean

setpoint_2 (*temperature=None, address=25*)
 Sets target temperature of furnace. Modbus address - 24

Parameters **temp** (*float, int*) – temperature in °C

Returns True if succesful, False if not

Return type Boolean

settings()

timer_end_type (*input=None, address=328*)

timer_resolution (*val=None, address=320*)

timer_status (*status=None, address=23*)

timer_type (*input=None, address=320*)

3.3 Module LCR

class `laboratory.drivers.lcr.LCR`

Bases: `object`

Driver for the E4980A Precision LCR Meter, 20 Hz to 2 MHz

| Attributes | message |
|----------------------|-------------------------------------|
| <code>maxtry</code> | max number to attempt command |
| <code>status</code> | whether the instrument is connected |
| <code>address</code> | port name |

| Methods | message |
|---------------------------|--|
| <code>connect</code> | attempt to connect to the LCR meter |
| <code>configure</code> | configures device for measurements |
| <code>write_freq</code> | transfers desired frequencies to the LCR meter |
| <code>trigger</code> | gets impedance for one specified frequency |
| <code>get_complexZ</code> | retrieves complex impedance from the device |
| <code>reset</code> | resets the device |

configure (*freq*)

Appropriately configures the LCR meter for measurements

display (*mode=None*)

Sets the LCR meter to display frequencies as a list

function (*mode='impedance'*)

Sets up the LCR meter for complex impedance measurements

get_complexZ ()

Collects complex impedance from the LCR meter

list_mode (*mode=None*)

Instructs LCR meter to take a single measurement per trigger

reset ()

Resets the LCR meter

shutdown ()

Resets the LCR meter and closes the serial port

trigger ()

Triggers the next measurement

write_freq (*freq*)

Writes the desired frequencies to the LCR meter

Parameters *freq* (*np.ndarray*) – array of frequencies

3.4 Module MFC

class `laboratory.drivers.mfc.AlicatController` (*port, address*)

Bases: `alicat.serial.FlowController`

Driver for an individual Mass Flow Controller.

Note: not called directly - access is from within :class:'~Drivers.MFC'

| Methods | message |
|--------------|--|
| get_massflow | gets massflow from controller |
| set_massflow | sets massflow on controller |
| get_pressure | gets pressure from controller |
| set_pressure | sets massflow on controller |
| get_temp | gets pressure from controller |
| get_vol_flow | gets volumetric flow from controller |
| get_setpoint | gets current set point from controller |
| reset | resets the device |

*see FlowController for more methods

get_all()

massflow (*value=None*)

Get or set the massflow of the appropriate flowmeter

Parameters *value* (*float*) – desired massflow value

pressure (*value=None*)

Get or set pressure of the appropriate flowmeter

Parameters *value* (*float*) – desired massflow value

reset()

sets the massflow to 0

setpoint()

Gets the current set point of the appropriate flowmeter

temperature()

Gets the temperature of the appropriate flowmeter :returns: gas temperature :rtype: float

volume_flow()

Gets the volumetric flow of the appropriate flowmeter

class `laboratory.drivers.mfc.MFC`

Bases: `laboratory.drivers.mfc.AlicatController`

Global driver for all Mass Flow Controllers

Note: see AlicatController for methods to control individual gases

| Attributes | message |
|------------|---|
| maxtry | max number to attempt command |
| status | whether the instrument is connected |
| co2 | controls for the Carbon Dioxide (CO2) controller |
| co_a | controls for the coarse Carbon Monoxide (CO) controller |
| co_b | controls for the fine Carbon Monoxide (CO) controller |
| h2 | controls for the Hydrogen (H2) controller |
| address | computer port address |

| Methods | message |
|--------------------------|---|
| <code>close_all</code> | closes all controllers |
| <code>connect</code> | attempt to connect to the LCR meter |
| <code>flush_all</code> | flushes data from the input/output buffer of all devices |
| <code>fugacity_co</code> | returns a ratio of CO ₂ /CO to achieve desired oxygen fugacity |
| <code>fugacity_h2</code> | returns a ratio of H ₂ /CO ₂ to achieve desired oxygen fugacity |
| <code>reset</code> | resets the device |

Example

```
>>> import Drivers
>>> mfc = Drivers.MFC()
>>> mfc.co2.get_massflow()
```

close_all()

Closes all flow controllers

flush_all()

Flushes the input? buffer of all flow controllers

fo2_buffer (*temp, buffer, pressure=1.01325*)

Calculates oxygen fugacity at a given temperature and fo2 buffer

| input options | type of fo2 buffer |
|---------------|---------------------------|
| 'QFM' | quartz-fayalite-magnetite |
| 'IW' | iron-wustite |
| 'WM' | wustite-magnetite |
| 'MH' | magnetite-hematite |
| 'QIF' | quartz-iron-fayalite |
| 'NNO' | nickel-nickel oxide |
| 'MMO' | molyb |
| 'CCO' | cobalt-cobalt oxide |

Parameters

- **temp** (*float, int*) – Temperature in u'°C'
- **buffer** (*str*) – buffer type (see table for input options)
- **pressure** (*float, int*) – pressure in bar (default: surface pressure)

Returns log10 oxygen fugacity

Return type float

fugacity_co (*fo2p, temp*)

Calculates the ratio CO₂/CO needed to maintain a constant oxygen fugacity at a given temperature.

Parameters

- **fo2p** (*float, int*) – desired oxygen fugacity (log Pa)
- **temp** (*float, int*) – temperature (u'°C)

Returns CO₂/CO ratio

Return type float

fugacity_h2 (*fo2p*, *temp*)

Calculates the ratio CO2/H2 needed to maintain a constant oxygen fugacity at a given temperature.

Parameters

- **fo2p** (*float*, *int*) – desired oxygen fugacity (log Pa)
- **temp** (*float*, *int*) – temperature (u°C)

Returns CO2/H2 ratio

Return type float

reset ()

Resets all connected flow controllers to 0 massflow

3.5 Module Motor

class `laboratory.drivers.motor.Motor` (*ports=None*)

Bases: `object`

Driver for the motor controlling the linear stage

| Attributes | message |
|-----------------------|--|
| <code>maxtry</code> | max number to attempt command |
| <code>status</code> | whether the instrument is connected |
| <code>home</code> | approximate xpos where <code>te1 == te2</code> |
| <code>max_xpos</code> | maximum x-position of the stage |
| <code>address</code> | computer port address |

| Methods | message |
|------------------------|---|
| <code>home</code> | move to the center of the stage |
| <code>connect</code> | attempt to connect to the LCR meter |
| <code>move</code> | moves the stage the desired amount in mm |
| <code>get_xpos</code> | get the absolute position of the stage |
| <code>set_xpos</code> | moves the stage the desired amount in steps |
| <code>get_speed</code> | get the current speed of the stage |
| <code>set_speed</code> | sets the movement speed of the stage |
| <code>reset</code> | resets the device |
| <code>test</code> | sends stage on a test run |

center ()

Moves stage to the absolute center

get_speed ()

Gets the current speed of the motor

Returns speed of motor

Return type float

get_xpos ()

Gets the current position of the stage

Returns x-position of stage

Return type str

home()
Moves furnace to the center of the stage ($x = 5000$)

move(*displacement*)
Moves the stage in the positive or negative direction
Parameters **displacement** (*float*, *int*) – positive or negative displacement [in mm]

reset()
Resets the stage position so that $xPos=0$

set_speed(*motor_speed*)
Sets the speed of the motor
Parameters **motor_speed** (*float*, *int*) – speed of the motor in mm/s

set_xpos(*xpos*)
Moves the linear stage to an absolute x position
Parameters **xpos** (*float*, *int*) – desired absolute position of stage in controller pulses

test()
Sends the motorised stage on a test run to ensure everything is working

3.6 Module Other

`laboratory.drivers.other.get_ports()`
Returns a list of available serial ports for connecting to the furnace and motor
Returns list of available ports
Return type list, str

`laboratory.drivers.other.load_instruments()`

`laboratory.drivers.other.reconnect(lab_obj)`

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

I

- `laboratory`, 4
- `laboratory.config`, 1
- `laboratory.drivers`, 16
 - `laboratory.drivers.daq`, 9
 - `laboratory.drivers.furnace`, 10
 - `laboratory.drivers.lcr`, 12
 - `laboratory.drivers.mfc`, 12
 - `laboratory.drivers.motor`, 15
 - `laboratory.drivers.other`, 16
- `laboratory.setup`, 1
- `laboratory.utils`, 8
 - `laboratory.utils.calibrate`, 5
 - `laboratory.utils.data`, 5
 - `laboratory.utils.loggers`, 7
 - `laboratory.utils.notifications`, 7
 - `laboratory.utils.plotting`, 7

A

AlicatController (class in *laboratory.drivers.mfc*),
12
append_data() (in module *laboratory.utils.data*), 6
append_data() (*laboratory.setup.Setup* method), 1
arrhenius() (*laboratory.utils.plotting.LabPlots*
method), 8

C

calculate_resistivity() (in module *labora-
tory.utils.plotting*), 8
center() (*laboratory.drivers.motor.Motor* method), 15
close_all() (*laboratory.drivers.mfc.MFC* method),
14
cole() (*laboratory.utils.plotting.LabPlots* method), 8
cond_fugacity() (*laboratory.utils.plotting.LabPlots*
method), 8
cond_time() (*laboratory.utils.plotting.LabPlots*
method), 8
configure() (*laboratory.drivers.daq.DAQ* method), 9
configure() (*laboratory.drivers.furnace.Furnace*
method), 10
configure() (*laboratory.drivers.lcr.LCR* method), 12
connect() (*laboratory.drivers.daq.DAQ* method), 9

D

DAQ (class in *laboratory.drivers.daq*), 9
Data (class in *laboratory.utils.data*), 5
data() (in module *laboratory.utils.loggers*), 7
debug (*laboratory.setup.Setup* attribute), 1
delayed_start (*laboratory.setup.Setup* attribute), 2
delayed_start (*labora-
tory.utils.notifications.Messages* attribute),
7
device_error (*labora-
tory.utils.notifications.Messages* attribute),
7
device_status() (*laboratory.setup.Setup* method),
2
display() (*laboratory.drivers.lcr.LCR* method), 12
dt_to_hours() (in module *laboratory.utils.plotting*),
8

F

filename (*laboratory.utils.data.Data* attribute), 5
find_center() (in module *labora-
tory.utils.calibrate*), 5
flush_all() (*laboratory.drivers.mfc.MFC* method),
14
flush_input() (*laboratory.drivers.furnace.Furnace*
method), 11
flush_output() (*labora-
tory.drivers.furnace.Furnace* method), 11
fo2_buffer() (*laboratory.drivers.mfc.MFC* method),
14
freq (*laboratory.utils.data.Data* attribute), 5
fugacity_co() (*laboratory.drivers.mfc.MFC*
method), 14
fugacity_h2() (*laboratory.drivers.mfc.MFC*
method), 14
function() (*laboratory.drivers.lcr.LCR* method), 12
Furnace (class in *laboratory.drivers.furnace*), 10
furnace_profile() (in module *labora-
tory.utils.calibrate*), 5

G

Gas (class in *laboratory.utils.data*), 5
gas() (*laboratory.utils.plotting.LabPlots* method), 8
get_all() (*laboratory.drivers.mfc.AlicatController*
method), 13
get_complexZ() (*laboratory.drivers.lcr.LCR*
method), 12
get_gas() (*laboratory.setup.Setup* method), 2
get_impedance() (*laboratory.setup.Setup* method),
2
get_ports() (in module *laboratory.drivers.other*), 16
get_Re_Im() (in module *laboratory.utils.plotting*), 8
get_speed() (*laboratory.drivers.motor.Motor*
method), 15
get_temp() (*laboratory.drivers.daq.DAQ* method), 9
get_temp() (*laboratory.setup.Setup* method), 2
get_thermopower() (*laboratory.drivers.daq.DAQ*
method), 10
get_thermopower() (*laboratory.setup.Setup*
method), 2

`get_voltage()` (*laboratory.drivers.daq.DAQ method*), 10
`get_xpos()` (*laboratory.drivers.motor.Motor method*), 15

H

`heating_rate()` (*laboratory.drivers.furnace.Furnace method*), 11
`home()` (*laboratory.drivers.motor.Motor method*), 15

I

`imp_diameter()` (*laboratory.utils.plotting.LabPlots method*), 8
`Impedance` (*class in laboratory.utils.data*), 6
`impedance_fit()` (*in module laboratory.utils.plotting*), 8
`index_temp()` (*in module laboratory.utils.plotting*), 8
`indicated()` (*laboratory.drivers.furnace.Furnace method*), 11

L

`lab()` (*in module laboratory.utils.loggers*), 7
`laboratory` (*module*), 4
`laboratory.config` (*module*), 1
`laboratory.drivers` (*module*), 16
`laboratory.drivers.daq` (*module*), 9
`laboratory.drivers.furnace` (*module*), 10
`laboratory.drivers.lcr` (*module*), 12
`laboratory.drivers.mfc` (*module*), 12
`laboratory.drivers.motor` (*module*), 15
`laboratory.drivers.other` (*module*), 16
`laboratory.setup` (*module*), 1
`laboratory.utils` (*module*), 8
`laboratory.utils.calibrate` (*module*), 5
`laboratory.utils.data` (*module*), 5
`laboratory.utils.loggers` (*module*), 7
`laboratory.utils.notifications` (*module*), 7
`laboratory.utils.plotting` (*module*), 7
`LabPlots` (*class in laboratory.utils.plotting*), 7
`LCR` (*class in laboratory.drivers.lcr*), 12
`leastsq_circle()` (*in module laboratory.utils.plotting*), 8
`list_mode()` (*laboratory.drivers.lcr.LCR method*), 12
`load_data()` (*in module laboratory.utils.data*), 6
`load_data()` (*laboratory.setup.Setup method*), 2
`load_frequencies()` (*laboratory.setup.Setup method*), 2
`load_instruments()` (*in module laboratory.drivers.other*), 16
`load_instruments()` (*laboratory.setup.Setup method*), 3

M

`massflow()` (*laboratory.drivers.mfc.AlicatController method*), 13
`Messages` (*class in laboratory.utils.notifications*), 7
`MFC` (*class in laboratory.drivers.mfc*), 13
`MFC_data` (*class in laboratory.utils.data*), 6
`Motor` (*class in laboratory.drivers.motor*), 15
`move()` (*laboratory.drivers.motor.Motor method*), 16

O

`other()` (*laboratory.drivers.furnace.Furnace method*), 11

P

`preflight_checklist()` (*laboratory.setup.Setup method*), 3
`pressure()` (*laboratory.drivers.mfc.AlicatController method*), 13

R

`read_errors()` (*laboratory.drivers.daq.DAQ method*), 10
`reconnect()` (*in module laboratory.drivers.other*), 16
`reconnect()` (*laboratory.setup.Setup method*), 3
`reset()` (*laboratory.drivers.daq.DAQ method*), 10
`reset()` (*laboratory.drivers.furnace.Furnace method*), 11
`reset()` (*laboratory.drivers.lcr.LCR method*), 12
`reset()` (*laboratory.drivers.mfc.AlicatController method*), 13
`reset()` (*laboratory.drivers.mfc.MFC method*), 15
`reset()` (*laboratory.drivers.motor.Motor method*), 16
`restart_from_backup()` (*laboratory.setup.Setup method*), 3
`run()` (*laboratory.setup.Setup method*), 3

S

`save_data()` (*laboratory.setup.Setup method*), 3
`save_obj()` (*in module laboratory.utils.data*), 6
`send_email()` (*in module laboratory.utils.notifications*), 7
`set_fugacity()` (*laboratory.setup.Setup method*), 3
`set_speed()` (*laboratory.drivers.motor.Motor method*), 16
`set_xpos()` (*laboratory.drivers.motor.Motor method*), 16
`setpoint()` (*laboratory.drivers.mfc.AlicatController method*), 13
`setpoint_1()` (*laboratory.drivers.furnace.Furnace method*), 11
`setpoint_2()` (*laboratory.drivers.furnace.Furnace method*), 11
`settings()` (*laboratory.drivers.furnace.Furnace method*), 11

Setup (*class in laboratory.setup*), 1
 shut_down() (*laboratory.setup.Setup method*), 4
 shutdown() (*laboratory.drivers.daq.DAQ method*), 10
 shutdown() (*laboratory.drivers.lcr.LCR method*), 12
 step_complete (*laboratory.utils.notifications.Messages attribute*), 7

T

Temp (*class in laboratory.utils.data*), 6
 temperature() (*laboratory.drivers.mfc.AlicatController method*), 13
 temperature() (*laboratory.utils.plotting.LabPlots method*), 8
 test() (*laboratory.drivers.motor.Motor method*), 16
 Thermo (*class in laboratory.utils.data*), 6
 timer_end_type() (*laboratory.drivers.furnace.Furnace method*), 11
 timer_resolution() (*laboratory.drivers.furnace.Furnace method*), 11
 timer_status() (*laboratory.drivers.furnace.Furnace method*), 11
 timer_type() (*laboratory.drivers.furnace.Furnace method*), 11
 toggle_switch() (*laboratory.drivers.daq.DAQ method*), 10
 trigger() (*laboratory.drivers.lcr.LCR method*), 12

V

voltage() (*laboratory.utils.plotting.LabPlots method*), 8
 volume_flow() (*laboratory.drivers.mfc.AlicatController method*), 13

W

write_freq() (*laboratory.drivers.lcr.LCR method*), 12